


DRUG-FREE PAIN RELIEF • RELAX OROFACIAL MUSCLES

The BNS-40

Portable Electrical Muscle Stimulator



The **BNS-40** is a battery operated, ultra-low frequency electrical muscle stimulator designed specifically for the dental patient. A personal, miniaturized version of the Myomonitor, it is ideally suited when prescribed for the following purposes:

- **Temporomandibular Disorders and associated pain**
- **Relieving symptoms associated with muscle spasm**
- **Relaxing orofacial muscles**
- **Increasing local blood circulation**
- **Increasing mandibular range of motion**

A time-tested, safe and effective means of drug-free relief for pain associated with muscle tension and a time-saving adjunct to TENS application in the dental office.

DRUG-FREE PAIN RELIEF • RELAX OROFACIAL MUSCLES

IN THE HOME . . .

- For drug-free treatment of the pain patient who is awaiting treatment or chooses not to undergo permanent treatment
- Save valuable office time and space by having the patient deprogram muscles at home prior to a bite registration appointment
- Effective relief of cervical pain associated with muscle tension

IN THE OFFICE . . .

- Save time — relax patients while the office model is in use for clinical procedures
- Relax spastic muscles after long procedures




INCREASE YOUR PROFIT OPPORTUNITIES

Whether the patient buys or rents the BNS-40, you will be assured of expanding your scope of services.

Use of the BNS-40 modality is often eligible for compensation under major medical plans. Its versatility in treating masticatory muscle pain associated with muscle spasm allows prescription over a wide range of circumstances.

“I find muscle hyperactivity to be a key element in most tension headache/TMD patients. I utilize the BNS-40 for these patients in the office and often prescribe it for home use to treat the muscular components. These patients usually experience dramatic relief, and acute episodes of the underlying conditions diminish. We have several units in our office and they are in use constantly.”

BARRY C. COOPER, DDS, NEW YORK, NY



The BNS-40

Pocket-sized, battery powered, and with unique true bilateral stimulation, it is the only instrument of its type specifically designed for personal use by the dental patient.

SPECIFICATIONS

Dimensions: 4.3" x 2.5" x 0.85"

Weight: 5 oz. with Battery & Patient Cables

Display: LED

Battery: 9-volt (alkaline) 250 hours

Unique Advantages:

- Pocket sized and lightweight
- Patient protection circuitry
- Optimized pulse width and frequency

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician or dentist.



MYOTRONICS
Helping Build the Perfect Bite

Phone (800) 426-0316 or (206) 243-4214

Fax: (206) 243-3625

www.myotronics.com • info@myotronics.com

